

SCISSORS

TUNGSTEN CARBIDE



STAINLESS STEEL

SUPERCUT

SCISSORS

Gorney Freeman F/L Straight

REF: GES-2010

TUNGSTEN
CARBIDE

Spread (open) shanks conform to the normal resting position of the hand thereby improving control and reducing fatigue.



Gorney Freeman F/L Curved

REF: GES-2020

TUNGSTEN
CARBIDE

Spread (open) shanks conform to the normal resting position of the hand thereby improving control and reducing fatigue.



Gorney F/L 5 5/8 In. Curved

REF: GES-2030

STAINLESS
STEEL

Original design with serrated stainless steel blades.



SCISSORS

Gorney F/L 5 5/8 In. Straight

REF: GES-2770

STAINLESS
STEEL

Original design with serrated stainless steel blades.



Gorney F/L 7 1/2 In. Straight

REF: GES-2040

STAINLESS
STEEL

Original design with serrated stainless steel blades.



Aston F/L Serrated 6 1/2 In. Curved

REF: GES-2050

TUNGSTEN
CARBIDE

Designed for undermining skin flaps. Double edge blades facilitate spread dissection.



SCISSORS

Aston F/L Serrated 8.0 In. Curved

GES-2060

TUNGSTEN
CARBIDE

Designed for undermining skin flaps. Double edge blades facilitate spread dissection.



Gorney Freeman F/L Curved

REF: GES-2070

TUNGSTEN
CARBIDE

Spread (open) shanks conform to the normal resting position of the hand thereby improving control and reducing fatigue.



Castanares F/L 6 1/4 In. Straight

GES-2080

TUNGSTEN
CARBIDE

Serrated with short sturdy blades for spreading soft tissue and cutting skin. Rounded tips prevent "buttonholing".



SCISSORS

PAR 4 1/2 In. Curved

GES-2090

TUNGSTEN
CARBIDE

Delicate with fine blades for precise dissection and cutting of delicate tissue.



PAR 5 3/4 In. Curved

GES-2100

TUNGSTEN
CARBIDE

Delicate with fine blades for precise dissection and cutting of delicate tissue.



PAR Serrated 4 1/2 In. Curved

GES-2110

TUNGSTEN
CARBIDE

Lower blade serrated to hold tissue or suture securely for clean, precise cut.



SCISSORS

PAR Serrated 5 3/4 In. Curved

GES-2120

TUNGSTEN
CARBIDE

Lower blade serrated to hold tissue or suture securely for clean, precise cut.



Stevens Tenotomy 4 1/2 In. Straight

GES-2130

TUNGSTEN
CARBIDE

Delicate blades with tapered tips that allow for precise cutting and dissection of fine tissues.



Stevens Tenotomy 4 1/2 In. Curved

GES-2140

TUNGSTEN
CARBIDE

Delicate blades with tapered tips that allow for precise cutting and dissection of fine tissues.



SCISSORS

Stevens Tenotomy 4 5/8 In. Curved

GES-2150

STAINLESS
STEEL

Excellent dissecting scissors for delicate tissue. Stainless steel blades.



Stevens Tenotomy 4 5/8 In. Straight

GES-2160

STAINLESS
STEEL

Excellent dissecting scissors for delicate tissue. Stainless steel blades.



Stevens Tenotomy 6 1/4 In. Curved

GES-2170

STAINLESS
STEEL

Excellent dissecting scissors for delicate tissue. Stainless steel blades.



SCISSORS

Devine 4 1/4 In. Curved

GES-2180

TUNGSTEN
CARBIDE

Fine dissecting scissors with sharp tips. Short, heavy blades minimize deflection when cutting or dissecting through tissue.



Strabismus 4 1/2 In. Curved

GES-2190

TUNGSTEN
CARBIDE

Delicate dissecting scissors with rounded tips. Useful for plastic and general dissecting surgery needs.



Aston F/L Serrated 6 1/2 In. Curved

REF: GES-2050

TUNGSTEN
CARBIDE

Designed for undermining skin flaps. Double edge blades facilitate spread dissection.



SCISSORS

Strabismus 4 1/2 In. Serrated Curved

GES-2210

TUNGSTEN
CARBIDE

One blade serrated to hold tissue or suture securely for clean precise cut.



Strabismus 4 1/2 In. Serrated Straight

GES-2220

TUNGSTEN
CARBIDE

One blade serrated to hold tissue or suture securely for clean precise cut.



Iris 4 1/2 In. Curved

GES-2230

TUNGSTEN
CARBIDE

Sharp tipped scissors provide superior distal cut of tissue and suture.



SCISSORS

Iris 4 1/2 In. Straight

GES-2240

TUNGSTEN
CARBIDE

Sharp tipped scissors provide superior distal cut of tissue and suture.



Iris 4 1/2 In. Serrated Curved

GES-2250

TUNGSTEN
CARBIDE

Lower blade serrated to hold tissue or suture for clean precise cut.



Iris 4 1/2 In. Serrated Straight

GES-2260

TUNGSTEN
CARBIDE

Lower blade serrated to hold tissue or suture for clean precise cut.



SCISSORS

Iris Supercut 4 1/2 In. Curved

GES-2270

TUNGSTEN
CARBIDE

Extremely sharp blades useful for cutting and dissection of tough tissues.



Iris Supercut 4 1/2 In. Straight

GES-2280

TUNGSTEN
CARBIDE

Extremely sharp blades useful for cutting and dissection of tough tissues.



Iris Swedish Pattern 4 1/4 In. Curved

GES-2290

STAINLESS
STEEL

General purpose Iris scissors with long blades and rounded shanks. Ideal for suture removal or delicate dissection.



SCISSORS

Iris Swedish Pattern 4 1/4 In. Straight

GES-2300

STAINLESS
STEEL

General purpose Iris scissors with long blades and rounded shanks. Ideal for suture removal or delicate dissection.



Iris American Pattern 4.0 In. Straight

GES-2310

STAINLESS
STEEL

General purpose Iris scissors ideal for suture removal or delicate dissection.



Iris American Pattern 4 1/2 In. Curved

GES-2320

STAINLESS
STEEL

General purpose Iris scissors ideal for suture removal or delicate dissection.



SCISSORS

Iris American Pattern 4 1/2 In. Straight

GES-2330

STAINLESS
STEEL



General purpose Iris scissors
ideal for suture removal or
delicate dissection.

McCollough Metzenbaum 5 3/4 In. Curved

GES-2340

TUNGSTEN
CARBIDE



Heavy pattern Metzenbaum.
Ideal for trimming tough
tissue.

Sealy Plastic 4 3/4 In. Curved

GES-2350

TUNGSTEN
CARBIDE



General purpose dissecting
scissors for precise cutting of
delicate tissues.

SCISSORS

Sealy Plastic 4 3/4 In. Straight

GES-2360

TUNGSTEN
CARBIDE

General purpose dissecting scissors for precise cutting of delicate tissues.



Septal Serrated 7 1/2 In. Angled

GES-2370

TUNGSTEN
CARBIDE

Designed for efficiency and flexibility in multiple surgical applications. Longer blades allow easy access.



Converse Angled 4.0 In.

GES-2380

STAINLESS
STEEL

Angled blades facilitate dissection of dense tissue.



SCISSORS

Nasal Tip Angled 3 3/4 In.

GES-2390

STAINLESS
STEEL



Very delicate angled blades facilitate spread dissection.

Joseph Nasal Curved 6 1/4 In.

GES-2400

STAINLESS
STEEL



Popular nasal scissors useful in dissecting and cutting dense tissues.

Fomon Dorsal Serrated 5 1/2 In. Angled

GES-2410

TUNGSTEN
CARBIDE



Angled shanks provide superior visibility.

SCISSORS

Metzenbaum 5 3/4 In. Curved

GES-2420

TUNGSTEN
CARBIDE

Popular dissecting scissors
with smooth rounded
scissors.



Metzenbaum 7.0 In. Curved

GES-2430

TUNGSTEN
CARBIDE

Popular dissecting scissors
with smooth rounded
scissors.



Metzenbaum Nelson 11.0 In. Curved

GES-2440

TUNGSTEN
CARBIDE

Popular dissecting scissors
with smooth rounded
scissors.



SCISSORS

Metzenbaum Serrated 5 3/4 In. Curved

GES-2450

TUNGSTEN
CARBIDE

Lower blade serrated to hold tissue or suture for a clean precise cut.



Metzenbaum Serrated 7.0 In. Curved

GES-2460

TUNGSTEN
CARBIDE

Lower blade serrated to hold tissue or suture for a clean precise cut.



Metzenbaum Serrated 9.0 In. Curved

GES-2470

TUNGSTEN
CARBIDE

Lower blade serrated to hold tissue or suture for a clean precise cut.



SCISSORS

Metzenbaum Supercut 5 3/4 In. Curved

GES-2480

TUNGSTEN
CARBIDE

Extremely sharp blades useful for cutting and dissecting tough tissues.



Metzenbaum Supercut 7.0 In. Curved

GES-2490

TUNGSTEN
CARBIDE

Extremely sharp blades useful for cutting and dissecting tough tissues.



Metzenbaum Supercut 9.0 In. Curved

GES-2500

TUNGSTEN
CARBIDE

Extremely sharp blades useful for cutting and dissecting tough tissues.



SCISSORS

Delicate Metzenbaum 5 3/4 In. Curved

GES-2510

TUNGSTEN
CARBIDE

Light and well balanced with delicate rounded tips.



Delicate Metzenbaum 7.0 In. Curved

GES-2520

TUNGSTEN
CARBIDE

Light and well balanced with delicate rounded tips.



Delicate Metzenbaum 9.0 In. Curved

GES-2530

TUNGSTEN
CARBIDE

Light and well balanced with delicate rounded tips.



SCISSORS

Delicate Metzenbaum Serrated 9.0 In. Curved

GES-2570

TUNGSTEN
CARBIDE

Lower blade serrated to hold tissue or suture for a clean precise cut.



Metzenbaum Sharp/Sharp 5 3/4 In. Curved

GES-2580

TUNGSTEN
CARBIDE

Sharp/Sharp tipped delicate scissors for precise cutting all the way to the tips.



Metzenbaum Sharp/Sharp 7.0 In. Curved

GES-2590

TUNGSTEN
CARBIDE

Sharp/Sharp tipped delicate scissors for precise cutting all the way to the tips.



SCISSORS

Mayo 5 1/2 In. Curved

GES-2600

TUNGSTEN
CARBIDE

For routine dissection and cutting of all types; sturdy, durable and able to cut the heaviest tissues.



Mayo 5 1/2 In. Straight

GES-2610

TUNGSTEN
CARBIDE

For routine dissection and cutting of all types; sturdy, durable and able to cut the heaviest tissues.



Mayo 6 3/4 In. Curved

GES-2620

TUNGSTEN
CARBIDE

For routine dissection and cutting of all types; sturdy, durable and able to cut the heaviest tissues.



SCISSORS

Mayo 6 3/4 In. Straight

GES-2630

TUNGSTEN
CARBIDE

For routine dissection and cutting of all types; sturdy, durable and able to cut the heaviest tissues.



Mayo Serrated 6 3/4 In. Curved

GES-2640

TUNGSTEN
CARBIDE

Lower blade serrated to hold tissue or suture for a clean precise cut.



Mayo Serrated 6 3/4 In. Straight

GES-2650

TUNGSTEN
CARBIDE

Lower blade serrated to hold tissue or suture for a clean precise cut.



SCISSORS

Mayo Supercut 6 3/4 In. Curved

GES-2660

TUNGSTEN
CARBIDE

Extremely sharp blades for cutting and dissection of tough tissues.



Mayo Supercut 6 3/4 In. Straight

GES-2670

TUNGSTEN
CARBIDE

Extremely sharp blades for cutting and dissection of tough tissues.



Operating Sharp/Blunt 5 1/2 In. Straight

GES-2680

TUNGSTEN
CARBIDE

General utility scissors. One tip sharp, one tip blunt for cutting gauze, suture, ect.



SCISSORS

Operating Sharp/Blunt 5 1/2 In. Curved

GES-2690

TUNGSTEN
CARBIDE

General utility scissors. One tip sharp, one tip blunt for cutting gauze, suture, ect.



Operating Sharp/Sharp 5 1/2 In. Straight

GES-2700

TUNGSTEN
CARBIDE

Sharp tipped, useful for cutting delicate tissues.



Operating Sharp/Sharp 5 1/2 In. Curved

GES-2710

TUNGSTEN
CARBIDE

Sharp tipped, useful for cutting delicate tissues.



SCISSORS

Operating Blunt/Blunt 5 1/2 In. Straight

GES-2720

TUNGSTEN
CARBIDE

General utility scissors useful for cutting gauze, sutures, ect.



Operating Blunt/Blunt 5 1/2 In. Curved

GES-2730

TUNGSTEN
CARBIDE

General utility scissors useful for cutting gauze, sutures, ect.



Littler 5.0 In. Curved

GES-2740

STAINLESS
STEEL

Fine dissecting scissors with hole in the distal part of the blade so suture can be threaded if needed.



SCISSORS

Spencer Stitch Removal 3 1/2 In. Straight

GES-2750

STAINLESS
STEEL

Delicate suture removal scissors ideal for office use. Crescent shaped tip for easy insertion under fine suture.



O'Brien Angled Stitch 3 3/4 In.

GES-2760

STAINLESS
STEEL

Suture removal scissors ideal for office use. Angled to the side blades with sharp points for ease of suture removal.

