

Aesthetic Laser

Laser Eye Protection

Order #. SERN-C



Order #. SESP



K.M.I.'s line of ocular laser shields, is the most complete line of sizes we know of. Sizes ranging from neonate up to various adult sizes. The all stainless steel shields are placed on the eye for protection during laser surgery. A high polished interior with a satin exterior finish to protect the operator from laser light "bounce back". The SESP's are available with metal pick-up arms, or rubber suction cups. (specify size when ordering). SERN-C's are available with suction cup pick-up only. The SERN-C's are in adult medium sizes only.