




MEDI CAPTURE

ANALOG TO DIGITAL IMAGES WITH USB TECHNOLOGY...

Order No: ES-MC19



MediCapture USB, Digital Medical Image Capture Device with USB Function

For use with nearly all video imaging devices (in Endoscopy, Dental, Ultrasound,...)

From analog medical-image-devices digital images can be easily captured and saved by the push of a button or by optional footswitch. Chosen images can be transferred through the integrated USB-Port to the USB-Flash-Drive-Stick and additionally can be saved. The Flash-Drive-Stick is suitable for the connection to all USB compatible devices like PC, Laptops, Mac or Webpad. Through the Flash-Drive-Stick integrated software, images can be saved and stored as JPEG, BMP, DICOM or PCX. In addition, it is possible to connect the MediCapture directly to your PC or Laptop via a supplied USB-Cable.

www.featherlift.com