

DISSECTORS

Agris Dingman

REF: ADD 14 inches (36cm)

STAINLESS
STEEL



Rounded, flattened tip enhances dissection. Can be modified with sharp tip if desired. Sold in pairs only.

MacCollum Dingman

REF: MDD 16 inches (41cm)

STAINLESS
STEEL




Designed to facilitate tissue dissection.

Sub-Glandular Paddle

REF: SBP

STAINLESS
STEEL



The sub-glandular breast paddle is a one piece instrument. The paddle portion of the device is made of surgical grade stainless steel and measures 25cm in length. The lateral edges of the paddle are formed in a manner to allow for easy dissection with minimal effort.