

LIGHT SOURCE

Halogen Light Source 150 Watts

REF: HL-2005-150

Multiple port turret that will accommodate Storz, Wolf, ACMI and Olympus fiber-optic light cables. A field replaceable bulb reduces down time. Specify 110 or 220 volts.



Fiber-optic Cables

FOC-SS-3518 / 3.5MMX180CM STORZ/STORZ FOC-SS-3523 / 3.5MMX230CM / STORZ/STORZ

FOC-AA-3518 / 3.5MMX180CM ACMI/ACMI FOC-AA-3523 / 3.5MMX230CM ACMI/ACMI

FOC-WW-3518 / 3.5MMX180CM WOLF/WOLF FOC-WW-3523 / 3.5MMX230CM WOLF/WOLF

FOC-SS-4818 / 4.8MMX180CM STORZ/STORZ FOC-SS-4823 / 4.8MMX230CM STORZ/STORZ

FOC-AA-4818 / 4.8MMX180CM ACMI/ACMI FOC-AA-4823 / 4.8MMX230CM ACMI/ACMI

FOC-WW-4818 / 4.8MMX180CM WOLF/WOLF FOC-WW-4823 / 4.8MMX230CM WOLF/WOLF

Headlight System

REF: H5000V

Lightweight design with a 3.5mm bifurcated light cable. Adjustable spot size. Custom comfort headband allows for extended wear without fatigue.

